



# Universal Bioink™ Kit Digital Light Processing (DLP) Bioprinting Protocol

This is a proposed DLP bioprinting procedure with recommended printing parameters for the **Universal Bioink™**. Please adjust this protocol according to your experimental needs. To maintain the sterility of the product, work under sterile conditions. For mixing and preparation of the **completely reconstituted Universal Bioink™** see the [Universal Bioink™ - Mixing Protocol \(Without Cells\)](#) or the [Universal Bioink™ - Mixing Protocol \(With Cells\)](#) at [voxcellbio.com](http://voxcellbio.com).

## Materials Required

- Completely reconstituted Universal Bioink™
- Tartrazine (100 mg in vial with magnetic stir bar)
- DLP bioprinter
- 365 or 405 nm UV light
- 1.25 mL Phosphate Buffered Saline (PBS)
- Magnetic stirring hotplate
- Pipette(s)
- Pipette tips

## DLP Bioprinting of Universal Bioink™

---

Note: The print parameters are suggestions and may vary depending on the specifications of the bioprinter being used. These instructions serve as guidelines. Please, adjust the settings as needed.

Adding **tartrazine** (also known as Acid Yellow 23) to the **completely reconstituted Universal Bioink™** prior to DLP printing is recommended. **Tartrazine** acts as a photoabsorber, significantly increasing the resolution of DLP bioprinting.

**If bioprinting with cells, follow the preparation of the completely reconstituted Universal Bioink™ according to the [Universal Bioink™ Kit - Mixing Protocol \(With Cells\)](#). For bioprinting without cells, follow the preparation of the completely reconstituted Universal Bioink™ according to the [Universal Bioink™ Kit - Mixing Protocol \(Without Cells\)](#).**

1. Prepare a 150 mM **tartrazine stock solution** by pipetting 1.25 mL of PBS into the **VoxCell tartrazine** vial. Using a stirring hotplate, gently stir (100-200 rpm) until a homogenous mixture is obtained. Store any unused **tartrazine stock solution** at room temperature in the absence of light.

Note: If printing in sterile conditions, use a syringe filter (0.22 µm pore size) to sterile filter the **tartrazine stock solution** prior to step 3.





- Using a stirring hotplate, heat the **completely reconstituted Universal Bioink™** (containing the **Bioink Photoinitiator: LAP**) to  $38 \pm 2$  °C while gently stirring (300-500 rpm) until a homogeneous mixture is obtained.
- Add 10 µL/mL of the **tartrazine stock solution** to the **completely reconstituted Universal Bioink™** to provide a final tartrazine concentration of 1.5 mM.
- Transfer the **completely reconstituted Universal Bioink™** containing **tartrazine** onto the DLP bioprinter platform using a pipette. The **completely reconstituted Universal Bioink™** containing **tartrazine** is now ready for bioprinting.

Note: Ensure that the build platform is heated to 37 °C during bioprinting to prevent thermal gelling.

- After the bioprinting procedure is complete, immediately hydrate the constructs with a buffer or media of choice.

## Printing Parameters

The recommended printing parameters are shown below. Please, adjust them to your experimental needs. Depending on the type of printhead being used, the first 2 layers may need higher exposure time and increased light intensity for build plate adhesion.

	Tartrazine Concentration (mM)	Exposure Time (s/layer)	Light Intensity (%)	Layer Height (µm)
Universal Bioink™	1.5	8	60	50

