



General Information

Storage

Universal Bioink™ in solid form should be stored at -20°C, in the absence of light, for a maximum of 6 months.

The **completely reconstituted Universal Bioink™** should be stored at 2 to 8 °C, in the absence of light, for a maximum of 10 days.

Always reseal all **Universal Bioink™** kit jars and vials after use.

Intended Use

For research purposes only.

Safety Information

Work in a ventilated area and use suitable personal protective equipment. For more information, please refer to the Safety Data Sheets.

Specifications

Product Name	Universal Bioink™ Kit
Appearance of Powder	White at 20°C
Reconstituted Volume	10 mL
Elastic Modulus*	80 ± 10 kPa
Sterilization Method	Filtration, Aseptic Processing
Sterility, USP <71>	Sterile

* Oscillatory testing via Anton Paar MCR 302 Rheometer at 23 °C

